



Because milk matters – bactotype Mastitis Test System for rapid, flexible and meaningful Mastitis monitoring and control

Carsten Schroeder, Daland C. Herrmann, Christine Gaunitz, Marco Labitzke, Oliver Sasse, Claudia Engemann

## MASTITIS


## Fast facts

- Mastitis is an inflammation of the udder
- Mastitis is common in dairy herds
- Mastitis is causing important economic losses worldwide
- Mastitis can't be eradicated -- but can be reduced to low levels by good management of dairy cows
- Of the several causes of Mastitis **only microbial infection is important**
- Bacteria, fungi, yeasts and possibly viruses can cause udder infection but the **main agents are bacteria**
- **Most common pathogens are *S. aureus*, *S. agalactiae*, *S. dysgalactiae*, *S. uberis* and *E. coli***
- Other Mastitis pathogens can cause occasional herd outbreaks




MASTITIS



## Mastitis pathogens



- 2010 study in collaboration with MBFG Wunsdorf and Boehringer Ingelheim
- Culture method used



## Mastitis pathogens



### Contamination!

Can be managed by improved hygiene or vaccination. Antibiotics not (always) necessary.

### Contagious

- *Staphylococcus aureus*
- *Streptococcus agalactiae*
- *Mycoplasma* spp.


### Infectious!

Requires strict disease control and antibiotic treatment and even culling.

### Environmental

- *Streptococcus uberis*
- *Escherichia coli* (*E. coli*)
- *Klebsiella pneumoniae*
- *Trueperella pyogenes* (previously *Arcanobacterium pyogenes*)
- *Staphylococcus epidermidis*
- *Pseudomonas aeruginosa*
- *Klebsiella oxytoca*
- *Corynebacterium bovis*
- *Enterobacter aerogenes*
- *Pasteurella* spp.
- *Brucella melitensis*
- *Proteus* spp.
- *Prototheca zopfii* (achlorophyllic algae)
- *Prototheca wickerhamii* (achlorophyllic algae)

## Detection methods



Bacterial culture	Real-time PCR (qPCR)
Broad range of detection	Fast, sensitive diagnostics
Semi-quantitative analysis	Quantitative analysis (C <sub>T</sub> value, curve)
Low costs for Agar plates	Multiplexing
Time and labour-intensive	Milk and pooled milk samples from disease control programs (e.g. <i>Brucellosis</i> ) can be used
Not suitable for sodium acide stabilised milk samples	Detection of other pathogens from milk (BVDV) possible
Some pathogens are difficult to culture (e.g. <i>Mycoplasma</i> )	Detects bacteria which are difficult to grow ( <i>Mycoplasma</i> )

## Rapid, flexible and meaningful Mastitis monitoring and control

Bovine Mastitis is the single most common and most costly disease of dairy cows.

Mastitis reduces milk production and raw milk quality, making it unsuitable for human consumption. Treating a cow herd for mastitis is costly, leads to milk withholding after treatment and may force producers to cull their dairy cows.

Reliably identifying pathogens that cause Mastitis is the cornerstone of targeted therapy strategies.

From surveillance to diagnostics, INDICAL's modular **bactotype Mastitis Test System** enables customers to implement different strategies for rapid and reliable Mastitis testing to meet their specific needs.

# BACTOTYPE MASTITIS SYSTEM

Manual extraction of DNA from milk  
IndiSpin Mastitis Kit  
(spin column technology)


OR

Automated extraction of DNA from milk  
IndiMag Mastitis Kit  
(magnetic bead technology)

# INDICAL


BIOSCIENCE

Tailored to your specific testing requirements



Contact our experts!

## DNA extraction from milk




Milk is a challenging sample type to extract DNA for PCR

Mastitis milk is even more challenging due to the wide range of sample conditions



## DNA extraction from milk



IndiField Mastitis Kit

For **manual DNA isolation**  
(spin-column/silica-membrane)



2

Optimized chemical lysis and  
silica-membrane extraction

- 1:30 h (24 samples)
- 6 h (96 samples)



Homogenize



Pathogen Lysis Microtubes  
+ TissuLyser II


Fast and efficient disruption  
of gram positive and gram  
negative bacteria.



IndiMag Mastitis Kit

For **automated DNA isolation**  
(magnetic beads)


2



Optimized chemical lysis and  
magnetic bead extraction

- 1:30 h (24 samples)
- 3 h (96 samples)

## DNA extraction from milk



**Pathogen  
Lysis Microtubes S**

+

**Buffer ML**

+



**400 µl milk sample**

**Homogenisation  
(TissueLyser II)**

**2 x 8 min (30 Hz)**

+

**Buffer MVL**



## bactotype Mastitis qPCR assays




For screening/detection:

- **Four** ready-to-use Multiplex real-time PCR Kits For analysis of quarter milk, bulk tank and pool samples
- PCR Kit-specific Master Mix (including IC system)
- PCR-Kit-specific Positive Control
- One Protocol for all bactotype Mastitis PCR Kits:
  - 5 min 95°C
  - 10 s 95°C
  - 30 s 57°C
 } 40x
- qPCR protocol ~ 1 hour
- Total time <3 hours (96 samples)

## Internal PCR Control

Our bactotype Internal PCR Control is based on a house-keeping gene and monitors correct DNA extraction and amplification

- Serves as indicator for elevated somatic cell counts:
- Low  $C_T$  of the bactotype Internal PCR Control = High somatic cell count



Bacterial culture (dilution series)



vs. Mastitis PCR (commercial competitor and INDICAL)

Mastitis Pathogen	Commercial competitor PP (magnetic beads extraction + detection)		INDICAL MagAttract magnetic bead extraction + bactotype Mastitis HP3 PCR Kit for detection	
	cfu/ml	C <sub>T</sub> value	cfu/ml	C <sub>T</sub> value
<i>Staphylococcus aureus</i>	900	36.91	180	37.21
<i>Streptococcus agalactiae</i>	122	36.42	122	36.21


- bactotype Mastitis Test System equally sensitive detecting *S. agalactiae*
- bactotype Mastitis Test System more sensitive detecting *S. aureus*

## Sensitivity for contagious Mastitis vs. commercial PCR PP

DNeasy Mastitis Mini Kit (extraction) +  
bactotype Mastitis HP3 PCR Kit (detection)

Results are equal (*S. agalactiae*) or better  
(*M. bovis*, *S. aureus*)

compared to commercial competitor PP



## Case study

### Rapid and reliable identification of *Mycoplasma bovis*

#### Procedure

In a field study, several animals in a herd of 600 dairy cows showed signs of clinical mastitis with unknown pathogenic cause.

60 milk pools (each of 10 animals) were generated, and DNA was extracted using the IndiMag Mastitis Kit. DNA samples were analyzed using the bactotype Mastitis HP2+ PCR Kit to detect *Mycoplasma*, *S. agalactiae* and *M. bovis* DNA. During the study, 12 animals succumbed to the infection.

The remaining 588 milk samples were individually tested by bacterial culture analysis by an audited and qualified milk analysis laboratory.


Note, bacterial culture can not distinguish between *Mycoplasma* and *M. bovis*.

#### Results:

At least 6.83 % of the herd tested *Mycoplasma* positive and at least 4.50 % tested *M. bovis* positive the by the bactotype Mastitis HP2+ PCR Kit (pool testing; Figure A)\*.

7.99% of the herd tested *Mycoplasma* positive when testing individual cows using bacterial culture (Figure B).

**Bacterial culture analysis confirm results obtained by the bactotype Mastitis HP2+ PCR Kit.**



\* To normalize results, one confirmed DNA positive milk pool was regarded as one *Mycoplasma* positive animal, even though more than one individual milk sample may have contributed to that positive pool.

## Testing pooled milk samples (*Mycoplasma bovis*) – positive animals contributing to the pool

SAMPLE	C <sub>T</sub> ( <i>M. bovis</i> )
<b>Pool 16</b>	<b>30.64</b>
16-01	-
16-02	35.55
16-03	37.65
16-04	-
16-05	-
16-06	26.37
16-07	-
16-08	-
16-09	-
16-10	-
<b>Pool 59</b>	<b>24.89</b>
59-01	-
59-02	-
59-03	27.15
59-04	21.96
59-05	-
59-06	-
59-07	-
59-08	-
59-09	-
59-10	27.03

Are pooled milk samples commonly used for monitoring programmes suitable for mastitis surveillance?

SAMPLE	DILUTION	C <sub>T</sub> ( <i>M. bovis</i> )
<b>Sample 59/ 04</b>	-	<b>22.99</b>
	1:10	26.74
	1:20	28.11
	1:30	29.13
	1:40	28.91
	<b>1:50</b>	<b>28.80</b>
<b>Sample 16/ 06</b>	-	<b>28.77</b>
	1:10	33.79
	1:20	37.16
	1:30	34.43
	1:40	34.58
	<b>1:50</b>	<b>35.14</b>
<b>Sample 16/ 02</b>	-	<b>34.19</b>
	1:10	-
	1:20	-
	1:30	-
	1:40	-
	1:50	-



Testing pooled milk samples (*S. aureus* and *S. agalactiae*)

SAMPLE	DILUTION	C <sub>T</sub> ( <i>Saureus</i> )
<b>Sample 217</b>	-	<b>25.18</b>
	1:10	28.23
	1:20	30.40
	1:30	30.00
	1:40	30.70
	<b>1:50</b>	<b>32.24</b>
<b>Sample UL-2</b>	-	<b>29.17</b>
	1:10	33.60
	1:20	33.66
	1:30	33.68
	1:40	33.71
	<b>1:50</b>	<b>33.95</b>
<b>Sample UL-1</b>	-	<b>33.13</b>
	1:10	35.28
	1:20	36.21
	1:30	36.04
	1:40	36.47
	<b>1:50</b>	<b>38.21</b>

SAMPLE	DILUTION	C <sub>T</sub> ( <i>Sagalactiae</i> )
<b>Sample 172</b>	-	<b>15.71</b>
	1:10	15.66
	1:20	15.99
	1:30	17.02
	1:40	16.66
	<b>1:50</b>	<b>16.87</b>
<b>Sample 177</b>	-	<b>26.41</b>
	1:10	30.46
	1:20	31.50
	1:30	32.14
	1:40	32.29
	<b>1:50</b>	<b>32.14</b>
<b>Sample 176</b>	-	<b>35.44</b>
	1:10	39.64
	1:20	36.88
	1:30	36.18
	1:40	36.59
	<b>1:50</b>	<b>36.25</b>

Positive results can be obtained from pooled milk samples (up to 50 animals) if at least one strongly positive animal (C<sub>T</sub> <32 in the individual milk sample) contributed to that pool. Milk pools containing only **one weak positive** sample (C<sub>T</sub> >34) may not be detected in a pool of 50 animals.

# bactotype

Mastitis Software

# bactotype Mastitis Software

PCR Resulting Program - clc

File Routine Calibration User Options Information Test

Lot Manager Sample plate layout Load Cycler Data Results

### Results


Plate ID: test3, File name: 170217HP2plate4 - Chart Data - Format 1.txt  
 Cycler name: hp301, User name: admin, Result Time: 02.11.2017 10:57:02

Sample ID	Assay	Info	Info Text	Analyt 1	Result 1	Analyt 2	Result 2	Analyt 3	Result 3	Analyt 4	Result 4	Valid	Valid Time	Valid To	To LIS Time	Comment
sample66	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[+++] 14,61	Internal Control	[+++] 26,53	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample67	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[+++] 20,48	Internal Control	[+] 30,08	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample68	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[+++] 22,54	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample69	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[+++] 24,20	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample70	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample71	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[+] 34,21	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample72	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[+] 33,12	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample73	HP2+	✓		Mycoplasma	[+++] 20,27	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[+++] 22,15	<input type="checkbox"/>	02.11.2017	admin		
sample74	HP2+	✓		Mycoplasma	[+] 30,21	S. agalactiae	[?] No ct	Internal Control	[+] 30,22	M. bovis	[+] 30,99	<input type="checkbox"/>	02.11.2017	admin		
sample75	HP2+	✓		Mycoplasma	[+] 36,18	S. agalactiae	[?] No ct	Internal Control	[+] 32,00	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample76	HP2+	✓		Mycoplasma	[+++] 17,36	S. agalactiae	[?] No ct	Internal Control	[+++] 21,00	M. bovis	[+++] 19,87	<input type="checkbox"/>	02.11.2017	admin		
sample77	HP2+	✓		Mycoplasma	[+++] 15,31	S. agalactiae	[?] No ct	Internal Control	[+] 32,00	M. bovis	[+++] 18,67	<input type="checkbox"/>	02.11.2017	admin		
sample78	HP2+	✓		Mycoplasma	[check Amplification]	S. agalactiae	[?] No ct	Internal Control	[+] 30,46	M. bovis	[+++] 12,80	<input type="checkbox"/>	02.11.2017	admin		
sample79	HP2+	✓		Mycoplasma	[+++] 14,43	S. agalactiae	[?] No ct	Internal Control	[+] 34,00	M. bovis	[+++] 16,62	<input type="checkbox"/>	02.11.2017	admin		
sample80	HP2+	✓		Mycoplasma	[+++] 26,72	S. agalactiae	[?] No ct	Internal Control	[+] 32,42	M. bovis	[+++] 28,45	<input type="checkbox"/>	02.11.2017	admin		
sample81	HP2+	✓		Mycoplasma	[+++] 12,07	S. agalactiae	[?] No ct	Internal Control	[+++] 17,44	M. bovis	[+++] 13,43	<input type="checkbox"/>	02.11.2017	admin		
sample82	HP2+	✓		Mycoplasma	[+++] 14,05	S. agalactiae	[?] No ct	Internal Control	[+] 22,00	M. bovis	[+++] 15,56	<input type="checkbox"/>	02.11.2017	admin		
sample83	HP2+	✓		Mycoplasma	[+++] 15,05	S. agalactiae	[?] No ct	Internal Control	[+++] 23,28	M. bovis	[+++] 16,65	<input type="checkbox"/>	02.11.2017	admin		
sample84	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample85	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
sample86	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
NC	HP2+	✓		Mycoplasma	[?] No ct	S. agalactiae	[?] No ct	Internal Control	[?] No ct	M. bovis	[?] No ct	<input type="checkbox"/>	02.11.2017	admin		
PC	HP2+	✓		Mycoplasma	[+++] 27,05	S. agalactiae	[+] 28,42	Internal Control	[+] 29,00	M. bovis	[+] 28,93	<input type="checkbox"/>	02.11.2017	admin		

admin 06.11.2017 14:17:26

- Results: pathogen, Ct value and rough quantification (+, ++, +++) will be given
- Information to run and possible faults (e.g. Internal Control not detected → extraction failed etc) will also be shown
- Print out and export functions will be implemented

# bactotype Mastitis Software



- Amplification curves can be viewed for each sample (pathogen and assay) to check for a specific highlighted fault or just to check the curves

## Summary

Bovine Mastitis is the single most common and most costly disease of dairy cows.

Mastitis reduces milk production and raw milk quality, making it unsuitable for human consumption. Treating a cow herd for mastitis is costly, leads to milk withholding after treatment and may force producers to cull their dairy cows.

Reliably identifying pathogens that cause Mastitis is the cornerstone of targeted therapy strategies.

From surveillance to diagnostics, INDICAL's modular **bactotype Mastitis Test System** enables customers to implement different strategies for rapid and reliable Mastitis testing to meet their specific needs.

## bactotype Mastitis Test System features and advantages



4x



- **Highly sensitive, fast and easy to perform**
- **2.5 to 3 hours runtime** from DNA extraction to results
- **Identifies difficult to grow bacteria**, like *Mycoplasma*
- **Optimized sample prep kits** for manual AND automated extraction of bacterial DNA from milk samples without compromising on quality
- **Four assays – one protocol:** combine and run up to four ready-to-use multiplex qPCR kits for targeted identification of Mastitis-causing pathogens
- Individual milk samples and pools of up to 50 milk samples can be tested
- **Advanced Internal PCR Control** for monitoring correct DNA extraction and amplification;
- **can also be used as indicator of somatic cell count in milk samples**
- **Data analysis made easy** with the bactotype Mastitis Software
- **Open platform** designed, tested and validated to integrate with most common qPCR thermal cyclers (AriaMx, Mx3005P, ABI 7500, Bio-Rad CFX96)
- **Supply certainty:** High production capacity and safety stocks
- **Exceptional technical support** by experienced veterinarians and molecular biologists

